
SOLDADORA BLISTER MANUAL FARMACEUTICA

PHARMACEUTICAL MANUAL BLISTER WELDER

Modelo: **SMC**

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COMPLEMENTOS

Útiles soldadores

WILSON GLOBAL INVEST offers an update on its tested and proven range of rugged manual welders for pharmaceutical sector of small-series heat and pressure welding of blister packs, with heat-sealable compounds like paper, aluminium or plastics, for to packing of capsules or other products.

OPERATION

- Use screws to adjust the height of the heatin plate with regard to the welding tool.
- Set the general switch to 1, lighting the green pilot light.
- Prefix the temperature on the temperature regulator.
- Wait 15 minutes until the welding plate reaches 120-150 °C and its ready to work when the green thermostat led lighs.
- Place the blister pack in the welding plate.
- Place capsules to be packed inside the blister.
- Place the welding tool under the welding plate.
- Lower the lever manually until it reaches its horizontal position.
- Wait 5-8 seconds (depending the aluminium film) and manually raise the lever to its vertical position.
- Remove the welding tool, remove the packaged pieces, and repeat the productive procces.

TECHNICAL SPECIFICATIONS

- Welding area 170 x 200 mm.
- Heating plate power 600 W
- Power supply 230 V F-N+TT
- Consumption 2,8A
- Temperature probe Type J. Removable by bayonet
- Pressure exerted 560 dan
- Manual Strength 4,5 dan on end of lever
- Recovery Via balanced anti-seizure spring
- Dimensions 570 mm. Height, 205 mm. Width, 370 mm. Depth 29 kg
- Weight

ACCESSORIES

Welding tools

WILSON GLOBAL INVEST, S.L.

Carretera Castellar 539 Pl. Ind. Can Petit 08227 Terrassa (España)

Tel. +34 93 93 664 57 32 - wilsongi@wilsongi.com - www.wilsongi.com

MAQUINARIA DE TERMOCONFORMADO
Y ENVASADO BLISTER O SKIN-PACK